
Patient chair CENTRIC excite


incl. short cover

plate 130mm

Electromotive height control, timeless classic design, high backrest, neck pad (head rest) adjustable via integrated holding magnets. Optional: Tiltable backrest, manual or electromotive seat adjustment, armrest pivotable backwards, foot rest and calf pad pivotable upwards, upper part of chair turnable 360° with detent at 45° and each 90°, durable elastic quality synthetic leather SKAI by Hornschuch in attractive colours appropriate for wiping disinfection.



850



HS DOMS GmbH

Langenschader Strasse 67 D-07318 Saalfeld

Fon: +49 (0) 36 71 / 54 50-0 Fax: +49 (0) 36 71 / 54 50-20

info@doms.eu



CENTRIC exciteType hp (chair shifting)



CENTRIC exciteCalf pad pivotable upwards (moving upwards synchronically with inclined backrest) and foot rest with foldable tread

 ϵ

Medical Device Directory / product class: 1 Certified according to Quality Management System DIN EN ISO 9001:2000

Technical Data:

DOMS TELECENTRIC excite incl. 2-instruments sliding table

Height: 1750 mm (w/o chart projector / LCD screen)

Length: 2570 mm Width: 1560 mm

1950 mm (2-instruments table moved into position 2)

Height of instruments table: 860mm

Patient chair Centric excite Minimal height of seat: 470 mm* Maximal height of seat: 670 mm* Lifting power: 1500 N

Connection voltage: 230V / 50Hz

Power: 1000 VA

Distributor:


Voltage supply (table instruments):

Low voltage: 6V / 12V

Non-heating connector: 230V / WS Protection (unit): IP 20 / Class: I

(*height of patient chair incl. seat upholstery Centric excite)

-dition 01/15



DOMS TELECENTRIC excite



